

## 2 Motor Massage and Heat System

### OWNER'S MANUAL and OPERATING INSTRUCTIONS

1. Before applying electrical power to your InSeat™ Massage System, place both UPPER and LOWER massage switches in the OFF position. Set the HEAT switch OFF, and the PULSE switch to ▼▼.
2. Plug cord from AC Transformer into mating receptacle extending from your InSeat™ product. Plug transformer body into wall outlet [120VAC 60Hz power].
3. Sit or lie in a comfortable position on your InSeat™ product.
4. Move UPPER massage switch to ON position. You will experience a steady massage action.
5. Push UPPER switch further to the PULSE position. Note the repetitive pulsing massage action.
6. Turn UPPER switch OFF, repeat steps 4 and 5 using the LOWER massage switch. The massage will now be felt on the opposite side of the lumbar area.
7. Now set both UPPER and LOWER switches to PULSE position. Both massage areas are now pulsing together. Move PULSE switch to ▲▼. Note the unique alternating action between the two message areas.
8. Turn both UPPER and LOWER switches OFF. Set HEAT switch ON. In a few minutes you will feel a comforting warmth in the lumbar area.

You now have completed all of the operational choices provided on your InSeat™ product. To indulge in the ultimate combination, select any of the above desired modes for massage and heat.

The massage and heat feature will automatically turn off after approximately 30 minutes of operation. To restart move any switch from OFF to ON. If all switches were on before timer turns off, move any switch to OFF, then back ON again.

*If you need assistance for any reason, please contact us:*

*Write to: Customer Service InSeat Solutions, Inc. Corporate Office  
11839 Smith Avenue, Santa Fe Springs, CA 90670*

*or call us Toll Free from anywhere in the United States at:*

**1-877-9-InSeat ext. 209**  
**(1-877-946-7328)**

Limited Manufacturer's Warranty

Patents Pending

©Copyright 2001 InSeat Solutions, Inc. Santa Fe Springs, CA

## IMPORTANT SAFETY INSTRUCTIONS

Read all instructions before using the InSeat™ System or the Relaxor® products.

### DANGER!

#### To reduce risk of electric shock.

1. Always unplug the product from the electrical outlet immediately after using and before cleaning.
2. Do not reach for a product that has fallen into water. Unplug immediately.
3. Do not use while bathing or in a shower.

### WARNING!

#### To reduce risk of burns, fire, electric shock, or injury to persons.

1. Your InSeat™ product should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
3. Never operate the InSeat™ product if it has a damaged cord or plug, or if it is not working properly. If it has been dropped or damaged, or dropped into water. Contact your dealer for instructions on returning it for examination and repair.
4. Close supervision is necessary when this appliance is used by, on, or near children or disabled persons.
5. Do not use outdoors.
6. Keep the cord away from heated surfaces.
7. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
8. If you are subject to any medical condition involving the use of medical device products, and/or if you are under a health professional or physician's care, please contact them before using this product.

**OPERATE  
YOUR  
SYSTEM  
ACCORDING  
TO THIS  
MANUAL**

**USE ONLY  
THE POWER  
DEVICE  
SUPPLIED  
WITH THIS  
PRODUCT**

### MAINTENANCE INSTRUCTIONS:

No user maintenance of this product is required. If the System does not function, be sure...

- (a) The TRANSFORMER or ADAPTOR is plugged into an electric outlet.
- (b) The TRANSFORMER or ADAPTOR cord is plugged into, and fully engaged in its receptacle.

If there is a malfunction of any kind, contact your dealer for instructions on returning it for examination and repair.

**SAVE THESE INSTRUCTIONS!**